

## Curved Electrospinning/Electrospraying needles



### Specification summery

Genvolt offers a variety of single, coaxial, triaxial, and bent needles featuring the industry-standard Luer-lock hub to cover all needs at your laboratory for research activities, for pre-industrial purposes. Genvolt has an increasingly prominent position in the market for complete range of needles used majorly in Laboratory & Research field.

Our angled needle system product is our 90-degree angled stainless stell reusable needles made of size 11 gauge.

### Needles chart

Gauge	Outside diameter (mm)	Inside diameter (mm)	Wall thickness (mm)	Needle Length (mm)
11	3	2.4	0.3	233

Highly functional needle system

We manufacture all needles using high quality materials.

Metal Chrome plate brass

Available in assorted sizes

## Schematic and dimensions

