


Side by Side Electrospinning/Electrospraying needles



Specification summery

Genvolt offers a variety of single, coaxial, triaxial, and bent needles featuring the industry-standard Luer-lock hub to cover all needs at your laboratory for research activities, for pre-industrial purposes. Genvolt has an increasingly prominent position in the market for complete range of needles used majorly in Laboratory & Research field.

Our side-by-side needle system product is our straight stainless steel reusable needles made of size 8 gauge.

Needles chart

Gauge	Outside	Inside	Wall	Triaxial	Arm
	diameter	diameter	thickness	needle	length
	(mm)	(mm)	(mm)	Length	(mm)
				(mm)	
8	4	3.2	0.4	81.85	22

Highly functional needle system

We manufacture all needles using high quality materials.

Metal Chrome plate brass

Available in assorted sizes



Schematic and dimensions

